



Contribution ID: 1072 Contribution code: THPM001

Type: **Poster Presentation**

Magnetic field control in the MedAustron synchrotron

Thursday 5 June 2025 15:30 (2 hours)

MedAustron, a synchrotron-based ion therapy centre in Austria, is focused on enhancing cancer treatment performance. A key improvement opportunity lies in the regulation of the main ring bending dipoles, which currently require time-consuming procedures to ensure reproducibility and reliability of the associated magnetic fields (B-Fields). Other therapy centres globally address this through a traditional B-train or similar systems to regulate on the B-Field and mitigate parasitic effects. In contrast to that, we propose a novel approach utilising a single Hall probe measurement inside a reference magnet, fused with a magnet model to provide real-time, high accuracy estimates of the integral B-Field for regulation. This technique, combined with a tailored controller, is evaluated under typical therapy cycling conditions. The system's performance is demonstrated through metrological analysis and beam property comparisons. Most importantly, the results show the possibility of significant improvements in treatment time reduction. Ultimately, the already achievable beam position accuracy, and spill structure in the treatment rooms, enable the start of commissioning in 2025.

Footnotes

Paper preparation format

LaTeX

Region represented

Europe

Funding Agency

Author: MARGREITER, Thomas (EBG MedAustron GmbH)

Co-authors: PLASSARD, Fabien (EBG MedAustron GmbH); PROKOPOVICH, Dale (EBG MedAustron GmbH); WOLF, Markus (EBG MedAustron GmbH); ADLER, Laurids (EBG MedAustron GmbH); PELLETIER, Sebastien (EBG MedAustron GmbH); GFALL, Immanuel (EBG MedAustron GmbH); BENEDIKT, Michael (European Organization for Nuclear Research)

Presenter: MARGREITER, Thomas (EBG MedAustron GmbH)

Session Classification: Thursday Poster Session

Track Classification: MC6: Beam Instrumentation and Controls, Feedback and Operational Aspects:

MC6.A28 Medical Applications